



Sardine freshness test

Slice fresh sardine in half and rinse each slice in tap water and Pursanova™ water for 3 minutes. Shake off the residual water adhered from rinsing and store at refrigerator for 3 days.



Tap Water Sardine Pursanova™
Water



Tap water rinsed sardine fillet started to show visible changes (fillet and skin) and odor occurred. 3 days later.

Pursanova™ water is maintaining its freshness and no odor detected. 3 days later.



Tap Water



Pursanova™ Water



Tap Water



Pursanova™ Water